


WWW.nsk-inc.com

NSK Offers

A Wide Choice of

VARIOS Ultrasonic Scaler Tips

For Clinical Applications.




Varios 2 / Varios tips

NSK tips fit Varios and Varios2 ultrasonic scaler handpieces as well as Satelec® ultrasonic scaler handpieces.

NSK make a number of tips which will fit on EMS® ultrasonic scaler handpieces.

The following is a guide to illustrate all NSK Varios tips, their applications and connectivity.



03 04

INTRODUCTION

This manual aims to help you choose ultrasonic tips and advises you on their use. For decades, ultrasonics have been used in odonto-stomatogy as well as in many surgical specialities. However, this technology is still underused in treatment in comparison with itsnumerous clinical applications.

NSK offers a range of "non-specialised" tips and specialised micro-tips.

Forincreased effectiveness, these micro-tips require a stand-alone generator with adjustable and reproducible power together with an independent irrigation system.

The system is broken down into 3 elements:

- The generator
- The handpiece
- The tip

The piezoelectric generator provides an accurate and reproducible power scale with 3 operating ranges. The different modes, P (perio), E (endo) and G (general) allow the power to be selected according to the tip's clinical indication.

The independent reservoirs enable irrigation by specific antiseptics at suitable temperatures.

NSK handpieces, whether optic or non-optic are lightweight and easy to handle. They provide improved tactile sensitivity in areas with no visual access. The smaller diameter, particularly in the optic handpiece, gives access to even the most posterior areas of the oral cavity.

There are many tips for dealing with various clinical indications.

They are designed to provide access and thus optimise the action of the working end fitting. Tip selection is guided by the treatment objective. (ref tables)



Application power range for NSK varios tips when used with NSK Varios2 utrasonic scaler



Easily interchangeable Solution bottles



Twin LED lignts

FUNDAMENTAL BASIS OF ULTRASONICS

The physics of ultrasonics on the same template as the piezosurgery handbook:

Ultrasonic equipment which combines irrigation with the movement of the tip in a single area produces 6 distinct effects:

Scanning: acts on the biofilm

Hammering: acts on mineralised deposits

Abrasion: acts on mineralised deposits and granulation tissue

Washing: acts on the biofilm and mineralised deposits

Micro-streaming: acts on the biofilm

Cavitation: acts on decontaminating the pocket

Depending on the treatment objectives, these effects are sought "more or less".


It is important to understand that the vibration movement, depending on its orientation on the work surface, has 2 clinical applications.

Scanning the biofilm when the tip is tangential to the root surface.

Hammering tartar when the tip is perpendicular to the debris to be removed.







SCALING

SCALING

MODEL G1

- · Removal of supra gingival tartar.
- · Removal of tartar narrow interdental space.
- Removal of marginal gingival calculus.





NSK Satelec® MODEL G1 EMS MODEL G1-E

Z217101 ORDER CODE Z291101

Single rooted / Multi rooted

NSK Satelec® MODEL G2 EMS® MODEL G2-E

MODEL G2

· Removal of supra gingival tartar.

23.4

Z291102

Cross Section Diagram

Z217102

Single rooted / Multi rooted

NSK Satelec® MODEL G3 EMS® MODEL G3-E

ORDER CODE Z291103

Single rooted / Multi rooted

* Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.

G3 MODEL

· Removal of supra gingival tartar.







Z217103

MODEL **G4**

- · Removal of supra gingival tartar.
- Removal of marginal gingival calculus.







Z217104 ORDER CODE Z291104

Single rooted / Multi rooted

MODEL G5

- · Removal of supra gingival tartar.
- · Removal of tartar narrow interdental space.
- · Removal of marginal gingival calculus.
- · Removal of gingival biofilm under natural tooth.







NSK Satelec® MODEL G5

Z217105

MODEL G6

- · Removal of supra gingival tartar.
- Removal of tartar narrow interdental space.
- · Removal of marginal gingival calculus.
- · Removal of gingival biofilm under natural tooth.







ORDER CODE <u>Z29</u>1106

Cross Section Diagram

Z217106





Refer to the power guide on pages 63-64.

Satelec® is registered trademarks of Acteon SA, France. EMS® is registered trademarks of Electro Medical Systems SA, Switzerland.

NSK Satelec® MODEL G6

EMS® MODEL G6-E

SCALING

MODEL G8

· Removal of supra gingival tartar.

Removal of tartar narrow interdental space.

· Removal of marginal gingival calculus.

22.7 Cross Section Diagram MODEL **G9**

· Removal of supra gingival tartar.

Removal of tartar narrow interdental space.

· Removal of marginal gingival calculus. · Removal of gingival biofilm under natural tooth.



MODEL G11

· Removal of supra gingival tartar.



NSK Satelec® MODEL G11

NSK Satelec® MODEL G11

Z217412

NSK Satelec® MODEL G11

Z217413

NSK Satelec® MODEL G8

Z217108

Single rooted / Multi rooted

NSK Satelec® MODEL **G9**

Z217109

Single rooted / Multi rooted

Z217411

Single rooted / Multi rooted

* Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.


MODEL G12

· Removal of supra gingival tartar.

MODEL G13

· Removal of supra gingival tartar.






Single rooted / Multi rooted



Single rooted / Multi rooted

* Satelec® is registered trademarks of Acteon SA, France.



PERIO

THE PHILOSOPHY OF ULTRASONIC TREATMENTS IN PERIODONTICS

Since the beginning of time, man has continually sought to push his boundaries. The desire for progress and the use of new technologies is particularly challenging for dentists who must now be able to integrate periodontics into their general practice. While conventional, non-surgical treatment is considered as the "Gold Standard", the fact remains that if access to a certain site is difficult, disinfection is not comprehensive and the quality of the curettage is still sometimes insufficient. These manual treatments are often tedious for the dentist and uncomfortable for the patient. Knowledge and understanding of ultrasonics and advances in the design of micro-tips are just some of the elements which mark the superiority of the ultrasonics treatment (Kocher, 1997 Gagnot et al in 1999 Gagnot, JPIO 2002) and push the boundaries of periodontal surgery. (Brent Scott et al 1999, Himeno 1994)(Ishikawa and Baehni, 2004; Heitz-Mayfield 2005)

We understand that the chances of reattachment and regeneration are related to the biocompatibility of tissues (cementum, gingival wall, periodontal ligament). Therefore the clinical challenge lies in the disruption of the biofilm and plaque retention elements (mineralised debris), while protecting tissues vital to healing (Goncalves and coll, 2008) (Fuzaka and Nishimura 1994. Smart et al 1990). Unlike other systems, the piezoelectric system provides a tip movement in the spindle of the handpiece. This movement controls the tip and allows the status of the surfaces to be respected by limiting iatrogenic effects. (Joki Kasthira et al 92) The preserved tissues will allow the formation of a new attachment. This control of ultrasonics and the movement of the tip allows us to tackle the idea of "tissue ecology", preserving the periodontal environment in order to permit long-term healing.

RECOMMENDATIONS FOR THE USE OF MICRO-TIPS IN PERIODONTICS

These usage recommendations aim to increase the effectiveness of the vibrating action and develop tactile sensitivity. These concepts are fundamental in optimising the effectiveness of the debridement.

ERGONOMICS

Adjustment in perio mode at low speed

Micro-irrigation adjustment.

Brush socket of the handpiece with three fingers: thumb on the NSK handpiece of the handpiece



Easing the weight of the cord Perioral pressure points



USING the tips

Light lateral pressure against the cemental surface

Work on the last two mm

Put in the salient line on the root surface (Viguier photos: file position 2mm from the tip)

Always remain in contact with the surface

Remember the concept of methodical, slow, limited scanning

Coronary access and lateral access of hard deposits.

These usage recommendations aim to increase the effectiveness of the vibrating action and

develop tactile sensitivity. These concepts are fundamental in optimising the effectiveness of

the debridement.

P5,6, P15~P18. writing by Drs Solenn HOURDIN(Ms) and Dominique GLEZ(Mr)

PRESENTATION OF THE RANGE OF MICRO-TIPS IN PERIODONTICS

The micro-tips used in the treatment of periodontitis (reference P) are different from the non-specialised tips (reference G) used on enamel for scaling.

The P micro-tips are used at low speeds in order to reduce iatrogenic effects on the cementum and to develop the dentist's tactile sensitivity.

■Scaling Tips ⇒ P.7-12

The profile of the tip's active area, located in the last 2mm, allows us to differentiate between them. A suitable profile makes it possible to work effectively at low speeds in deep pockets without injected anaesthetic.

There are three main families of micro-tip: Curette micro-tips, round micro-tips and diamond micro-tips (ref Gilles Gagnot book). These tips are used carefully and methodically. (ref tables)

Curette micro-tips (P10, P26R, P26L) are used in the initial stages of treatment. Round (P1, P20, P21R, P21L, P25R, P25L) and composite (V-P10) micro-tips are appropriate for the maintenance stage and diamond micro-tips (P1D, P2D, P3D) are kept for the treatment of abscesses, persistent inflammatory and furcation areas. ■Perio Tips ➡ P.19-28

The table presents the different indications for each of the tips and their usage speed. This guide gives details of the tips in combination with clinical videos in order to improve the learning curve. The equipment offered by NSK is suited to the implementation of current cleaning techniques by fully meeting the specifications.

					OCa	ling tip				
TIP No.	G1	G2	G3	G4	G5	G6	G8	G11	G12	G13
Mode selection	G	G	G	G	G	G	G	G	G	G
Single rooted		0	0	0	0		0	0	0	0
Multirooted		0	0	0	0	•	0		0	0
Removal supra gingival tartar	0	0	0	0	0	0	0	0	0	0
Removal tartar interdental narrow space					0	•	0			
Removal marginal gingival calculus				0	0		0			
Removal tartar in deep gingival										
Removal gingival biofilm					0					
Removal on Biofilm furcation area										
Removal of Implant Biofilm										
Removal on junction biofilm prosthetic tooth										
Removal on biofilm dental erosion										
Incising of furcations										
Incising of abscess										
Peeling internal wall of pocket										

	Micro tip round			Micro tip Curette			Micro tip diamonded			
TIP No.	P1	P20	P21R/L	P25R/L	V-P10	P10	P11R/L	P26R/L	P1D	P2D/3D
Mode selection	Р	Р	Р	Р	Р	Р	Р	Р	Р	Р
Single rooted	0	0	0	0	0	0	0	0	0	0
Multirooted		0			0	0		0	0	
Removal supra gingival tartar						0				
Removal tartar interdental narrow space	0					0	0	0	0	0
Removal marginal gingival calculus						0				
Removal tartar in deep gingival							0	0	0	0
Removal gingival biofilm	0	0				0				
Removal on Biofilm furcation area				0			0	0		
Removal of Implant Biofilm					0					
Removal on junction biofilm prosthetic tooth					0					
Removal on biofilm dental erosion					0					
Incising of furcations							0	0		0
Incising of abscess									0	0
Peeling internal wall of pocket									0	0

BIBLIOGRAPHY

Gonçalves et Coll. Root Cementum May modulate gene expression During Periodontal Regeneration: A preliminary study in Humans. J Periodontol., 2008. 323-331. Gilles Gagnot et Coll. Les ultrasons en odontologie. Applications Thérapeutiques. Memento Editions CdP

KOCHER T., RÜHLING A., MONSEN H., PLAGMANN HC. -Effectiveness of subgingival instrumentation with power-driven instruments in the hands of experienced and inexperienced operators. J Clin Periodontol 1997; 24: 498-504.

GAGNOT G., MICHEL JF., DARCEL J., CATHELINEAU G. - Action de nouveaux inserts ultrasoniques sur les dômes des espaces interradiculaires. Observation au MEB. J Parodontol Implant Orale 1999: 19: 411-419

GAGNOT G., MORA F., POBLETTE M., VACHEY E., MICHEL JF., CATHELINEAU G. -A comparative study of manual and ultrasonic instrumentation on root surfaces. J Periodontol 2002; 73: 10-1245.

BRENTSCOTT J., STEED-VEILANDS AM., YUKNA RA. -Efficacité accrue du détartrage des espaces inter-radiculaires par l'utilisation d'inserts ultrasonores diamantés. Parodontie et Dentisterie Restauratrice 1999; 19: 354-361.

HIMENO H., TSUTSUMIDA K., SATO K., HONDA M., MIYAHARA A. -A study on a new and practical ultrasonic root planning method (URM). The journal of the hokkaido Dental Association 1991; 46: 122-133.

FUKAZAWA E., NISHIMURA K. -Superficial cemental curettage: its efficacy in promoting improved cellular attachment on human root surfaces previously damaged by peridontitis, J Periodontol 1994: 65: 168-176.

SMART GJ., WILSON M., DAVIES EH., KIESER JB. -The assessment of ultrasonic root surface debridement by determination of residual endotoxin levels. J Clinical Perio 1990; 17: 174-178.

JOTIKASTHIRA NE., LIE T., LEKNES KN. -Comparative in vitro studies of sonic, ultrasonic and reciprocating scaling instrument. J Clin Periodontal 1992; 19 (8): 560-9.

Heitz-Mayfield LJ.

How effective is surgical therapy compared with nonsurgical debridement? Periodontal 2000, 2005:37:72-87

Ishikawa I. Baehni P.

Nonsurgical periodontal therapy--where do we stand now?

Periodontol 2000. 2004;36:9-13.



MODEL P1

- · Removal of tartar narrow interdental space.
- Removal of gingival biofilm.







ORDER CODE Z217404

MODEL P1D

- · Removal of tartar narrow interdental space.
- · Removal of tartar in deep gingival.
- · Debridement of abscess.
- · Peeling wall internal pocket.







NSK Satelec* MODEL P1D

Z217405

· Diamond coated



MODEL P2D

MODEL P3D

- · Removal of tartar narrow interdental space. · Removal of tartar in deep gingival. · Debridement of furcations.
- · Debridement of abscess. · Peeling wall internal pocket.

27.0














· Diamond coated · Left curved type





Cross Section Diagram

Z217403







* Satelec® is registered trademarks of Acteon SA, France.



MODEL P10

- · Removal of supra gingival tartar.
- Removal of tartar narrow interdental space.
- Removal of marginal gingival calculus.
- Removal of gingival biofilm.







NSK Satelec® MODEL P10

Z217410

Sharp edge

21



MODEL P11R

MODEL P11L

- Removal of tartar narrow interdental space. Removal of tartar in deep gingival. Removal of gingival biofilm.
- Removal of Biofilm furcation area. Debridement of furcations.






33.0









· Sharp edge · Left curved type

NSK Satelec® MODEL P11L



Z217409

MODEL P20

· Removal of gingival biofilm.

NSK Satelec® MODEL P20

Refer to the power guide on pages 63-64.

P20

Cross Section Diagram

Z217420

Single rooted / Multi rooted

ORDER CODE Z291420









































Z217422

Single rooted / Multi rooted

^{*} Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.





- · Removal of gingival biofilm.
- Removal of Biofilm furcation area.







Z217424 ORDER CODE Z291424

· Right curved type






MODEL P25L



Single rooted / Multi rooted


P11R / P11L / P21R / P20 / P21L (From above)



MODEL P26R

MODEL P26L

- · Removal of tartar narrow interdental space. · Removal of tartar in deep gingival. · Removal of gingival biofilm.
- Removal of Biofilm furcation area. Debridement of furcations.








NSK Satelec® MODEL P26L

· Sharp edge · Left curved type



Z217427

P25R / P25L / P26R / P26L (From above)




NSK Satelec® MODEL P26R

Refer to the power guide on pages 63-64.



Z217426





[·] Sharp edge · Right curved type

^{*} Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.

PERIO (Maintenance Tip)

MODEL V-10

MODEL V-P10

For the maintenance of prostheses and implants. Produced from special material to allow safe debridement without causing damage to metals and implants.

MODEL V-P11R MODEL V-P11L

V-P11R/L left and right are double bent tips. This tip is extra long to easily access difficult to reach areas and to avoid interference with neighbouring teeth or contours. Its shape was designed to best fit different tooth morphologies, making this a versatile tip for a wide range of uses from anterior teeth to premolars.















NSK Satelec® MODEL V10

Z217041

· Includes E-Tip replacement wrench * Plastic Tip is not included

· Pack of 3

* V10 holder is not included

MODEL V-P10

Y900184

MODEL V-P11R ORDER CODE Y1002165

· Pack of 3 · Right curved type

* V10 holder is not included

· Pack of 3 · Left curved type * V10 holder is not included

MODEL V-P11L



MODEL V-P12

It serves a wide range of applications, including the maintenance of implants and other prosthetic devices. Its shape offers easy access to all treatment areas.

V-P26R


V-P12 is a longer version of the popular V-P10 tip. V-P26R/L is a double-bent tips. The wide angles offer access to areas that cannot be reached with V-P11R/L tips or other conventional instruments. These tips are ideally suited to access the distal surfaces of molars and corner areas as, the end of the tip fits tooth morphology very closely.


MODEL V-P26L



18.1









MODEL V-P12

· Pack of 3 * V10 holder is not included



MODEL V-P26R ORDER CODE Y1002168

· Pack of 3 · Right curved type * V10 holder is not included

· Pack of 3 · Left curved type * V10 holder is not included

MODEL V-P26L ORDER CODE Y1002169

V10 holder: 1 pc.

 V-P11R, V-P11L, V-P12, V-P26R, V-P26L: 1 pc. each · Includes E-Tip replacement wrench

Perio-Control Tips

Starter Kit



Y1002166

ORDER CODE

^{*} Satelec® is registered trademarks of Acteon SA, France.


Refer to the power guide on pages 63-64.

PERIO (Irrigation and Cleaning)

For effective cleaning following debridement or scaling, *No scaling effect,

MODEL P41 MODEL P40

· Safe to use regardless of the pocket condition due to the marking guideline of 3 mm from the position of 6 mm.



Tip diameter ø 0.6 mm

Tip diameter ø 0.8 mm





MODEL Perio Kit A

MODEL Perio Kit B





NSK Satelec® MODEL Perio Kit A

Tip holder: 1 pc.

• P1D、P2D、P3D、P10: 1 pc. each

Y900144

MODEL Perio Kit B

Y900145

• P20、P25R、P25L: 1 pc. each • Tip holder: 1 pc.



ENDODONTICS

ENDODONTICS



MODEL **E4**

· Application for root canal cleaning and removal of foreign particles.



Tip diameter ø 0.6 mm

NSK Satelec® MODEL E4

MODEL **E7**

· Water not required

· Removal of filling and foreign material.

24.1

E8 MODEL

· Removal of filling and foreign material.



· Water not required



Cross Section Diagram

Z217307

MODEL **E4D**

· Application for root canal cleaning and root canal enlargement.



· Diamond coated

MODEL E7D

· Enlargement of root canal wall.




· Diamond coated

NSK Satelec® MODEL E7D

MODEL E8D

· Enlargement of root canal wall.



NSK Satelec* MODEL E8D

with diamond

28.4

· Diamond coated

Z217317



* Satelec® is registered trademarks of Acteon SA, France.



Cross Section Diagram

Z217428

* Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.

Refer to the power guide on pages 63-64.

Z217318



MODEL **E5**

· Lateral condensation.

MODEL **E6**

· Lateral condensation.

21.5

ORDER CODE Z217305 NSK Satelec* MODEL E5

· Water not required





· Water not required





MODEL E15D

· For locating root canals and removal of calcifications.






NSK Satelec® MODEL E15D

Z217319


· Diamond coated



* Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.







MODEL Endo Kit E11 For anterior teeth MODEL Endo Kit E12 For anterior teeth





NSK Satelec®

MODEL Endo Kit E11

Y900146

• E11 : 1 pc. 120 angle holder for anterior teeth

· E Tip replacement wrench : 1 pc.

· U files 33 mm #15, #20, #25, #30, #35 : Pack of 6 each size

NSK Satelec®

MODEL Endo Kit E12


Y900147

• E11 : 1 pc. 120 angle holder for anterior teeth

· E Tip replacement wrench : 1 pc.

• U files 33 mm #15, #20, #25, #30, #35 : Pack of 6 each size

		#15	#20	#25	#30	#35
E11	Varios2	22 mm	24.5 mm	22 mm	22 mm	22 _{mm}
E 11	Varios	22-23 mm	22-23 mm	23-25 mm	23-25 mm	24-26 mm
E12	Varios2	22 mm	24.5 mm	24.5 mm	24.5 mm	22 _{mm}
EIZ	Varios	21-25 mm	21-25 mm	21-25 mm	21-25 mm	21-25 _{mm}



* Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.

Refer to the power guide on pages 63-64.

* Satelec® is registered trademarks of Acteon SA, France.

The microscopic diamond coated tip allows effective removal of dentin at the apex

2 types of angled tips (90°, 70°) are available for anterior teeth and another 2 types of curved tips (right and left) are available for molar teeth.

MODEL E30RD

MODEL E30LD

· Diamond coated tip allows efficient cavity preparation.



· Diamond coated · For posterior teeth (right angled) · Diamond coated · For posterior teeth (left angled)





MODEL E31D

MODEL E32D

MODEL Retro Kit

· Diamond coated tip allows efficient cavity preparation.



Section

Z217210











Y900149

· Diamond coated · For anterior teeth (90°)

NSK Satelec® MODEL E32D

• E30RD, E30LD, E31D, E32D : 1 pc. each • Tip holder : 1 pc.



* Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.

* Satelec® is registered trademarks of Acteon SA, France.

NSK Satelec* MODEL E31D

· Diamond coated · For anterior and posterior teeth (70°)

EXTRACTION




Periodontal membrane removal Tip

MODEL G95

• Tooth extraction is possible without excessively cutting the alveolar bone around the remaining root, retaining the morphology of the alveolar bone.







- Use of irrigation enhances performance.
- When used with Varios 970 choice of irrigation.
- · Tip can be used on NSK ultrasonic units.

NSK MODEL G95

ORDER CODE Z217705

• TiN coating • 0.7mm thick • A double edge, which is sharp in three directions

Root extraction (upper left No. 2 tooth removal)





Impacted tooth (lower left No. 8 tooth removal)





Be careful not to apply excessive force, as in cases of using elevators.

Refer to the power guide on pages 63-64.

* Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.



RESTORATIVE



MODEL G61D

MODEL G62D

· Use the tip with the spherical diamond at the top. Bleeding occurs rarely.















· Diamond coated

NSK Satelec® MODEL G61D

Z217601

NSK Satelec* MODEL G62D

ORDER CODE Z217602

· Diamond coated



* Satelec® is registered trademarks of Acteon SA, France. Refer to the power guide on pages 63-64.

MODEL G65D

· For molar, interdental proximal caries. The hemispherical top of the tip does not allow cutting of sound tooth.





ORDER CODE Z217605 NSK Satelec* MODEL G65D

Diamond coated



MODEL G66D

· For anterior proximal caries. Only decayed dentin can be removed as per for molars.







NSK Satelec® MODEL G66D

ORDER CODE Z217606

· Diamond coated



MODEL G67D

· Diamond coated

· For molar, interdental proximal caries. The hemispherical top of the tip does not allow cutting of sound tooth.

MODEL G68D

 For anterior proximal caries. Only decayed dentin can be removed as per for molars.

















ORDER CODE Z217608

· Diamond coated

NSK Satelec® MODEL G68D





* Satelec® is registered trademarks of Acteon SA, France. Refer to the power guide on pages 63-64.



Condensation / Loosening / Plugging

MODEL G21

· For cementing a small prosthesis such as an inlay for anterior teeth.









· For cementing a large prosthesis such as an onlay or a crown for molar teeth.





MODEL Tip Cap



MODEL Tip Cap

Y140751

• For G21 / G22 • Pack of 10 • Autoclavable [135]c

MODEL G26

Post removal. (ø1.3 mm thickness)



· Amalgam condensation.









NSK Satelec® MODEL **G21**

• For anterior teeth • Water not required • 3 caps supplied

ORDER CODE Z217010

NSK Satelec® MODEL G22

• For anterior teeth • Water not required • 3 caps supplied

Z217020

Z217206


ORDER CODE Z217208

· Water not required





* Satelec® is registered trademarks of Acteon SA, France. Refer to the power guide on pages 63-64.



V-TIP SYSTEM

NSK Varios V-Tip holder system is innovative, multi-functional and economical.

The V-Tip holder is a separate component to the actual cutting tips and can be easily and inexpensively replaced,

NSK V-Tip holder system offers the flexibility and expandability to meet the wide requirements for multi-purpose clinical applications.



V-TIP System Line Up



3 Features

Safety

- · Minimal damage to soft tissue
- Less risk of bleeding
- Cutting with ultrasonic instruments is minimally invasive

Economical

- · One tip Replaceable inserts
- Inserts are inexpensive

Ease of Use

- · Small tip head guarantees excellent visibility
- · Small head enables great access to distal sides of molar teeth and buccal area
- One tip meets the wide requirements of many clinical applications




MODEL V-30 Restorative Pages 55~58 · For Minimal Intervention · For Finishing / Trimming / Polishing · For Caries Removal Z217030 • Includes E Tip replacement wrench *Tip is optional Endodontics Pages 53~54 **V30** 120° Angle type · For Root Canal Cleaning · For Root Canal Enlargement

51



For Root Canal Cleaning



V-TIP SYSTEM (Endodontics)

MODEL V-S51

Y900170

• For root canal cleaning 23 mm • Pack of 3

MODEL V-S52

Y900171

· For root canal cleaning 26 mm · Pack of 3 MODEL V-S53 Y900172

• For root canal cleaning 32 mm • Pack of 3





For Root Canal Enlargement

Refer to the power guide on pages 63-64.



• Diamond coated • For root canal enlargement 17 mm • Pack of 3

Y900173




Y900174

• Diamond coated • For root canal enlargement 23 mm • Pack of 3



Y900175 · Diamond coated · For root canal enlargement 26 mm · Pack of 3



• Diamond coated • For root canal enlargement 32 mm • Pack of 3



*In order to reduce the incidence of ledge formation, there is no diamond coating on the end of the tip. (1mm)



Removal of broken files

* Satelec® is registered trademarks of Acteon SA, France.

The advantage of ultrasonic cutting is that tip movement is minute allowing only decayed dentine to be removed, while sound dentine is preserved. Cervical decay can be treated with minimum bleeding enabling quicker recovery.



For Minimal Intervention



Y900118

· Diamond coated · Pack of 3

MODEL V-G71



Y900119

· Diamond coated · Pack of 3

MODEL V-G72



Y900120

· Diamond coated · Pack of 3





Cervical decay removal Root decay removal

For Minimal Intervention

MODEL V-G75



Y900121

· Diamond coated · Pack of 3

MODEL V-G76



Y900122

· Diamond coated · Pack of 3

Contact surface of tooth decay removal




Root decay removal

* Satelec® is registered trademarks of Acteon SA, France.

The NSK Varios V-Tip System has been developed to carry out "Minimally Invasive" dental procedures. As an altenative to rotary cutting, the Varios V-Tip offers precise and effective ultrasonic technologycutting.




For Finishing / Trimming / Polishing



· Includes E Tip replacement wrench *Tip is optional



For Caries Removal







MODEL Assorted Tip Kit
(For Minimal Intervention)

MODEL Assorted Tip Kit
(For Finishing / Trimming / Polishing)

MODEL Assorted Tip Kit
(For Caries Removal)







MODEL Assorted Tip Kit ORDER CODE

Y900158

MODEL Assorted Tip Kit ORDER CODE Y900168

MODEL Assorted Tip Kit ORDER CODE Y900164

. V-G70, V-G71, V-G72, V-G75, V-G76: 1 pc. each

• V-G77, V-G78, V-G79 : 1 pc. each

· V-S1, V-S3, V-S33 : 1 pc. each











Compact Multipurpose Control Unit at Your Fingertips



Periodontic, Endodontic, Scaling, MI applications and Implant maintenance

- Potable ultra compact body
 - Simple and intuitive operation control panel

· Powered by iPiezo engine

Specifications

· Large water adjustment and power knob for easy operation

: 100-240V AC 50/60 Hz

: W80 x D115 x H32 (mm)

(with out knob and cord)

: 11 W (G mode)

· Unit weight : 0.43 kg (without attachment)

- Holders are available for various installation
- Handpiece is washable in the thermodisinfector and autoclavable up to 135°C



Varios 370 Complete Sets

Optic Varios 370 Optic Complete Sets

MODEL VA370 LUX (120V)

Y1001956

MODEL VA370 LUX (230V)

Y1001343

Non-Optic Varios 370 Non-Optic Complete Sets

MODEL VA370 (120V)

Y1001958

MODEL VA370 (230V)

Y1001340

· Control unit with handpiece cord · Control unit holder · AC adaptor · Handpiece

· Handpiece holder · Water filter set · Water supply connector (120V only)

• Foot control (FC-71) • 3 Scaling tips (G4, G6, G8) • Tip wrench with torque limiter

Multifunction Portable Unit for Extreme Versatility

- All-in-Stylish compact all-in one unit for
- Stable vibration and power by newly developed iPiezo engine
- Convenient large capacity bottles of 400ml each or different solutions with independent pumps
- Excellent illumination of the treatment area is guaranteed by the twin LEDs guided through twin glass rods
- · Handpiece can be repeatedly autoclaved up to 135°C and processed in a thermo disinfector
- Extensive tip range for a wide range of applications
- Easily Interchangeable Solution Bottles
- Simple Unit Design for Clean Maintenance

Power Supply : AC 120V or 230V 50/60 Hz

: 11 W (G mode)

: W160 x D270 x H190 (mm)



Varios 970 Complete Sets

Varios 970 Optic Complete Sets

VA970 LUX (120V)

Y1001166

VA970 LUX (230V)

Y1001168

Non-Optic Varios 970 Non-Optic Complete Sets

MODEL VA970 (120V)

Y1001174

VA970 (230V)

Y1001175

- · Control unit · Handpiece · Handpiece cord · 3 Scaling tips (G4, G6, G8)
- Foot control (FC-70) 2 x Irrigation bottles 3 x Tip wrench with torque limiter
- · AC power cord · Sterilization case

· Unit weight : 2.1 kg (Except attachment) · Bottle volume : 400 mL (x 2 bottles)

POWER GUIDE

TIP NO.	Varios 2 970/370	Varios 750/560/350	Satelec® P-Max XS P5 XS	EMS® Piezon Master 600 (miniMaster)
G1	G 7	G 3	15	P6
G2	G 7	G Max	15	P6
G3	G 7	G 3	15	P3
G4	G 5	G 3	15	P3
G5	G 6	G 3	15	-
G6	G 8	G 1	15	P6
G8	G 7	G 3	15	-
G9	G 5	G 1	15	-
G11	G 6	G 3	15	-
G12	G 7	G 3	15	-
G13	G 6	G 3	15	-
	PERIC)/MAINTEN	IANCE	
P1	P 6	P 2	5	-
P1D	P 10	P 2	5	-
P2D	P 8	P 4	5	-
P3D	P 8	P 4	5	-
P10	P 10	P Max	5	-

	NSK		Satelec®	EMS®
TIP NO.	Varios2 970/370	Varios 750/560/350	P-Max XS P5 XS	Piezon Master 600 (miniMaster)
P11R	P 10	Max	5	-
P11L	P 10	Max	5	-
P20	P 8	P 4	5	P6
P21R	P 8	P 4	5	P6
P21L	P 8	P 4	5	P6
P25R	P 10	P 4	5	P 6
P25L	P 10	P 4	5	P6
P26R	P 6	P 3	5	-
P26L	P 6	P 3	5	-
P40	P 3	P 1	5	-
P41	P 3	P 1	5	-
V-P10	P 6	P 3	5	-
V-P11R	P 6	P 3	5	-
V-P11L	P 6	P 3	5	-
V-P12	P 6	P 3	5	-
V-P26R	P 6	P 3	5	-
V-P26L	P 6	P 3	5	-
				_

NSK Satelec® EMS®										
TIP NO.	Varios2	Varios 2 Varios		EMS® Piezon Master 600 (miniMaster)						
ENDODONTICS										
E4	6	3 3	5	-						
E4D	3 7	3 3	5	-						
E 5	1 10	3 3	5	-						
E6	1 10	3 4	5	-						
E7	6	⊜ Max	5	-						
E7D	B 6	⊜ Max	5	-						
E8	3 2	3 4	5	-						
E8D	3 2	3 4	5	-						
E15D	3 6	3 3	5	-						
E11	3 4	3 3	5	-						
E12	3	3 3	5	-						
E30RD	B 6	B 3	5	-						
E30LD	3 6	3 3	5	-						
E31D	3 6	3 3	5	-						
E32D	3 6	3 3	5	-						
	Е	XTRACTIO	N							
G95	G 10	_	_	_						

TIP NO.		NSK		EMS®						
111 110.	Varios2 970/370	Varios 750/560/350	P-Max XS P5 XS	Piezon Master 600 (miniMaster)						
G61D	G 2	G 2	15	-						
G62D	G 2	G 3	15	-						
G65D	G 2	G 2	15	-						
G66D	G 2	G 2	15	-						
G67D	G 2	G 3	15	-						
G68D	G 2	G 3	15	-						
	PI	ROSTHETIC	CS							
G21	G 4	G 2	20	-						
G22	G 4	G 2	20	-						
G26	G 10	G Max	20	-						
G28	G 10	G 2	20	-						
	V-TIP	(ENDODON	ITICS)							
V-S50	3 5	3 3	10	-						
V-S50D	3 5	3 3	10	-						
V-S51	3 3	1	10	-						
V-S51D	3 3	1	10	-						
V-S52	3 5	3 2	10	_						

	N:	SK	Satelec®	EMS®
TIP NO.	Varios2 970/370	<i>Varios</i> 750/560/350	P-Max XS P5 XS	Piezon Master (miniMaste
V-S52D	B 5	3 2	15	-
V-S53	B 5	3 2	15	_
V-S53D	B 5	3 2	15	-
	V-TIP(RESTORAT	IVE)	
V-G70	G 6	© 3	15	-
V-G71	G 6	© 3	15	-
V-G72	G 6	© 3	15	-
V-G75	G 6	G 1	15	-
V-G76	G 6	G 3	15	-
V-G77	G 6	G 3	15	-
V-G78	G 6	G 3	15	-
V-G79	G 6	G 3	15	-
V-S1	G 6	G 3	15	-
V-S3	G 4	© 3	15	-
V-S33	G 6	G 3	15	_

: Figures in the power guide indicate the appropriate power

range and maximum power level for each tip. Do not exceed the maximum power level to avoid possible tip breakage. Applicable power range for NSK Varios tips when used with NSK Varios and Varios2 ultrasonic scaler.

Varios2

-	911												
ı	POW	ER R	ange								gine*	⁄ari:	os2
	Gó				0	0	€	4	Ġ	G 0	8	9	0
	Εō		0 6	9 0 G	6 0 8	9 0							
	Ρó	006	0 60	9890									T
	0%	10%	20%	30%	40%	50%	. 6	50%	70%	80%	90	%	100%

Varios



<Caution>

- Always consult the operating manual and attached sheets before use.
- Always use with correct power mode and power range for each individual tip.
- Always use with sufficient irrigation to prevent excessive heat generation during treatment.
- Always start out with lower power setting and then adjust to higher setting if necessary.
- Always make sure that tips are securely tightened with the appropriate wrench and the power is set to and does not exceed the appropriate level before use.


*Incorrect operation may cause tip breakage even if tip is used within the recommended power range. Consult the operating manual before use.

Satelec* & P-Max are registered trademarks of Acteon SA, France. EMS*, Piezon & miniMaster are registered trademarks of Electro Medical Systems SA, Switzerland.

63

The Varios Tip Card - Easy to use guide to assess tip wear

Ultrasonic scaler's tips wear through repeated use. Worn tips swing differently to new tips and do not perform as effectively. To check the condition of the tip, use the Varios Tip Card.



How to use the Varios Tip Card

- 1. Place the neck of the tip in the cut out.
- 2. Check wear of the tip.
- 3. See the green, yellow and red line to check wear of the tip "See below what each color means."

At NSK we recommend to replace a tip when the tip meets the yellow line (wear of 1mm) to guarantee safe and effective use.





Tip Replacement Wrench with Torque Limiter

For tips G1- G68D, G95, P1- P41, E4, E4D, E15D, E30RD, E30LD, E31D, E32D and V-tip replaceable insets.

*Please tighten the tip until it clicks.



E Tip Replacement Wrench

For V10 (V-tip holder), tips E5, E6, E7, E7D, E8, E8D, E11, E12 and the inserts of V30 (V-tip holder).



NAKANISHI INC. www.nsk-inc.com 700 Shimohinata, Kanuma, Tochigi 322-8666, Japan

NSK Europe GmbH

Elly-Beinhorn-Strasse 8 65760 Eschborn, Germany

NSK Oceania Pty Ltd

Unit 22, 198-222 Young St. Waterloo, Sydney, NSW 2017, Australia

NSK United Kingdom Ltd

Office 5, Gateway 1000, Arlington Business Park Whittle Way, Stevenage, SG1 2FP, UK

NSK Asia

1 Maritime Square, #09-33 HarbourFront Centre 099253. Singapore

NSK America Corp

1800 Global Parkway, Hoffman Estates IL 60192. USA

NSK Middle East

Room 6EA701, 7th Floor, East Wing No.6 Dubai Airport Free Zone, PO Box 54316, Dubai, UAE

NSK America Latina Ltda

Dr. Placido Gomes, 610-Sala 104, A. Garibaldi Joinville SC 89202-050. Brazil

Specifications are subject to change without notice.

The Varios Tip Videos

Please watch the videos about the Varios tip range at the below websites to experience the precision and quality that is synonymous with NSK.

Ultrasonic scaler Varios series



■ Series description [Varios 970] ■ Maintenance Varios tip



■Tip variations
■Tip description

Doctors Comments



▶ ▶ ▶ NSK Web site

INTERNATIONAL	http://www.nsk-dental.com	SPAIN	http://www.nsk-spain.es
GERMANY	http://www.nsk-europe.de	RUSSIA	http://www.nsk-russia.ru
UNITED KINGDOM	http://www.nsk-uk.com	USA / CA	http://www.nskdental.com
FRANCE	http://www.nsk.fr	Latin America & Caribbean	http://www.latin-america.nsk-dental.com

Specifications are subject to c