



SURGICAL HANDPIECES

X-SG20L
X-SG25L
X-SG93L/X-SG93
X-SG65L

accurate & efficient



The NSK Surgical Handpiece range offers Clinicians the comfort and minute dexterity for all precise surgical applications. Excellent visibility direct to the operational field is greatly enhanced via the integrated cellular optic illumination. Smooth and precise torque with a high load capability is delivered even at very low speeds.



Powerful Partners®

Ti-Max X Surgical Handpieces

- Titanium Body with Scratch Resistant DURACOAT
- Cellular Glass Optics (X-SG20L / X-SG25L / X-SG93L / X-SG65L)

X-SG20L

20:1 Reduction



External and internal cooling available



Optic MODEL X-SG20L ORDER CODE C1003

- Push Button Chuck • Sealed Head Design
- External and internal cooling (Kirschner and Meyer)
- Max. Speed 2,000 min⁻¹

X-SG25L

1:1 Direct Drive



Optic MODEL X-SG25L ORDER CODE C1011

- Clean Head System • Push Button Chuck
- External cooling • Max. Speed 40,000 min⁻¹

X-SG93L / X-SG93

1:3 Increasing



Triple Spray



Optic MODEL X-SG93L ORDER CODE C1004

Non-Optic MODEL X-SG93 ORDER CODE C1007

- Clean Head System • NSK "UP2" Push Button Chuck
- External cooling • Max. Speed 120,000 min⁻¹

X-SG65L

Straight Handpiece

1:1 Direct Drive



Optic MODEL X-SG65L ORDER CODE H1009

- For HP burs (ø2.35) • For CA burs (ø2.35)*
- *Need to insert the supplied bur stopper
- Clean Head System • Max. Speed 40,000 min⁻¹



Washable in the thermos disinfectant



Autoclavable up to max 135°C

NAKANISHI INC.
www.nsk-inc.com

700 Shimohinata Kanuma-shi
Tochigi 322-8666,
Japan
TEL: +81 (0)289-64-3380
FAX: +81 (0)289-62-5636

NSK Europe GmbH EC REP
www.nsk-europe.de

Elly-Beinhorn-Strasse 8
65760 Eschborn,
Germany
TEL: +49 (0)6196 77606-0
FAX: +49 (0)6196 77606-29

NSK United Kingdom Ltd
www.nsk-uk.com

Office 5, Gateway1000,
Arlington Business Park, Whittle Way,
Stevenage, SG1 2FP, UK
TEL: +44 (0)1438 310670
FAX: +44 (0)1438 310671

NSK Oceania Pty Ltd
www.nsk-inc.com

Unit 22, 198-222 Young St.
Waterloo, Sydney,
NSW 2017, Australia
TEL: +61 (02) 8306 3000
FAX: +61 (02) 8306 3030

NSK America Corp
www.nsk-inc.com

700 Cooper Court
Schaumburg, IL 60173,
USA
TEL: +1 847 843 7864
FAX: +1 847 843 7822

NSK Middle East
www.nsk-inc.com

Room 6EA-701, 7th Floor, East Wing No.6
Dubai Airport Free Zone,
PO Box 54316 Dubai, UAE
TEL: +971 (0)4-204-5507
FAX: +971 (0)4-204-5508